

---

# **tomolog-cli Documentation**

***Release 0.1***

**Argonne National Laboratory**

**Nov 03, 2022**



## CONTENTS



<b>1 Installation</b>	<b>3</b>
<b>2 Usage</b>	<b>15</b>
<b>3 API reference</b>	<b>17</b>
<b>4 Credits</b>	<b>21</b>
<b>Bibliography</b>	<b>23</b>
<b>Python Module Index</b>	<b>25</b>
<b>Index</b>	<b>27</b>



**tomolog-cli** is a command-line-interface for publishing tomography experiment data and meta data in a stack of Google slides

### L\_var\_plum\_038

- File name: L\_var\_plum\_038.h5
- Beamline: 2-BM Micro-tomography
- Scan date: 2022-11-02T17:27:34-0500
- Exposure time: 0.70000 s
- Camera pixel size: 3.45  $\mu\text{m}$
- Lens magnification: 7.5
- Projection pixel size: 0.46  $\mu\text{m}$
- Angle step: 0.120 °
- Number of angles: 1501 (0.00 - 180.00)
- Projection size: 3232 x 2426
- Scan energy: 25.51 keV
- Sample Y: 35.79 mm
- Propagation dist.: 30.00 mm
- Load Raw: 11.33325 V
- Load: -11620.56799 N



Frame from the IP camera in the hutch



```
tomocupy recon --file-name /data/2022-10/Stock/L_var_plum_038.h5 --remove-stripe-method fw --reconstruction-type full --rotation-axis-auto auto --find-center-end-row 1500
```



---

**CHAPTER  
ONE**

---

## **INSTALLATION**

First, you must have [Conda](#) installed and create a dedicated conda environment:

```
(base)$ conda create -n tomolog python=3.9
```

and:

```
(base)$ conda activate tomolog
(tomolog)$
```

then install all [requirements](#) with:

```
(tomolog)$ conda install -c conda-forge python=3.9 dropbox google-api-python-client
          ↵matplotlib dxchange dxfile python-dotenv opencv matplotlib-scalebar
```

install meta

```
(tomolog)$ git clone https://github.com/xray-imaging/meta.git
(tomolog)$ cd meta
(tomolog)$ python setup.py install
```

and install tomolog

```
(tomolog)$ git clone https://github.com/xray-imaging/tomolog-cli.git
(tomolog)$ cd tomolog
(tomolog)$ python setup.py install
```

### **1.1 Requirements**

Please install all the packages listed in [requirements](#) file.


**tomolog** also requires access tokens from dropbox and google services.

### **1.1.1 Dropbox**

Go to [dropbox developer site](#) , login using your google credentials and select “Create an App”:



Take the App key and App secret from the Settings tab:



The screenshot shows the 'Creating a Dropbox app' section of the Dropbox developer console. It includes three numbered steps:

- ① Configure app settings**: Name your app and choose initial settings.
- ② Select access scopes**: Choose the access scopes, or specific permissions, that your app needs to interact with Dropbox. We recommend starting small and adding more permissions later if you need them. [Get started](#)
- ③ Add branding**: Give your users important information about your Dropbox app. Should comply with the Dropbox developer branding guide. [Get started](#)

Below this, there are sections for Status (Development selected), Development teams (0 / 1), Development users (Only you), Permission type (Scoped App), OAuth 2 Redirect URIs, Allow public clients (Implicit Grant & PKCE), Generated access token, Chooser / Saver / Embedder domains, Webhooks, and Extensions.

an copy them in a file in your home directory called:

```
$ ~/.tomologenv
```

```
APP_KEY=....  
APP_SECRET=....
```


Set the following permissions:


The screenshot shows the "Creating a Dropbox app" section of the developer console. It includes three steps: 1. Configure app settings (Name your app and choose initial settings), 2. Select access scopes (Choose the access scopes, or specific permissions, that your app needs to interact with Dropbox. We recommend starting small and adding more permissions later if you need them. [Get started](#)), and 3. Add branding (Give your users important information about your Dropbox app. Should comply with the Dropbox developer branding guide. [Get started](#)). Below this, there are sections for Status (Development), Development teams (0 / 1, with buttons for Enable additional teams and Unlink all teams), Development users (Only you, with a button for Enable additional users), Permission type (Scoped App), OAuth 2 Redirect URIs (https:// (http allowed for localhost), Add), Allow public clients (Implicit Grant & PKCE) (Allow dropdown set to Allow), Generated access token (Generate button), Chooser / Saver / Embedder domains (example.com, Add button), and Webhooks Webhook URIs (https://, Add button). At the bottom, there is an Extensions section with an Extensions URI (https://).

### 1.1.2 Google


Next step is to authorize tomolog to create slides on the ...@gmail.com gmail account.

Open a web browser and login as ...@gmail.com then go to [google developer site](#) and press “Select a project” to create a new project







Once the project is selected go to the project dashboard:



and select API & Services / Credentials



select it and then go to Create Credentials / Service account



where you enter the service account name, i.e. tomolog

The screenshot shows the Google Cloud IAM & Admin Service Accounts creation interface. On the left, a sidebar lists various IAM components: IAM, Identity & Organization, Policy Troubleshooter, Policy Analyzer, Organization Policies, Service Accounts (which is selected and highlighted in blue), Workload Identity Federation, Labels, Tags, Settings, Privacy & Security, Identity-Aware Proxy, Roles, Audit Logs, and Manage Resources. At the top, the title is "Create service account". The main area is divided into three steps:

- Service account details**: Fields include "Service account name" (placeholder: "Display name for this service account"), "Service account ID" (placeholder: "Email address: <id>@user7bmb-365420.iam.gserviceaccount.com"), "Service account description" (placeholder: "Describe what this service account will do"), and a "CREATE AND CONTINUE" button.
- Grant this service account access to project (optional)**
- Grant users access to this service account (optional)**

At the bottom of the main area are "DONE" and "CANCEL" buttons.

Grant this service account access to project: Owner

Grant users access to this service account:

```
service account user role: Google Account email: ....@gmail.com
service account admin role: Google Account email: ....@gmail.com
```

Once the service account is selected you need to create the authorization key:

The screenshot shows the Google Cloud IAM & Admin Service Accounts Keys tab for the "tomolog" service account. The sidebar on the left is identical to the previous screenshot. The main area has tabs: DETAILS, PERMISSIONS, KEYS (which is selected and highlighted in blue), METRICS, and LOGS. The KEYS tab displays a section titled "Keys" with a warning message: "Service account keys could pose a security risk if compromised. We recommend you avoid down can learn more about the best way to authenticate service accounts on Google Cloud [here](#)." Below this is a note: "Add a new key pair or upload a public key certificate from an existing key pair." and links to "organization policies" and "Learn more about setting organization policies for service accounts". A "ADD KEY" button is at the bottom.

## Create private key for "test"

Downloads a file that contains the private key. Store the file securely because this key can't be recovered if lost.


- Key type
- JSON  
Recommended
- P12  
For backward compatibility with code using the P12 format

CANCEL      [CREATE](#)


open the json file and copy the authorization token in a file called:

```
$ ~/tokens/google_token.json
```

Go back to the project dashboard




and select “Go to the API overview”



and select “Enable API and Services”



search and add the “google slide api”



1 result



### Google Slides API

Google Enterprise API [?](#)

Reads and writes Google Slides presentations.

Now tomolog will be able to create slides in the ... @gmail.com user account.



---

## CHAPTER TWO

---

## USAGE

To publish experiment log information to a google page:

```
$ tomolog run --file-name /local/data/2022-03/Peters/B4_Pb_03_c_10keV_892.h5 --  
→presentation-url https://docs.google.com/presentation/d/  
→128c8JYiJ5EjbQhAtegYYetwDUVZILQjZ5fUIoWuR_aI/edit#slide=id.p
```

For help:

```
$ tomolog run -h  
usage: tomolog run [-h] [--file-name PATH] [--PV-prefix PV_PREFIX] [--beamline {None,2-  
→bm,7-bm,32-id}] [--idx IDX] [--idy IDY] [--idz IDZ] [--max MAX] [--min MIN]  
[--presentation-url PRESENTATION_URL] [--rec-type {recgpu,rec}] [--  
→config FILE] [--config-update] [--double-fov] [--logs-home FILE] [--token-home FILE]  
[--verbose]  
  
optional arguments:  
-h, --help            show this help message and exit  
--file-name PATH      Name of the hdf file (default: .)  
--PV-prefix PV_PREFIX          PV prefix for camera (default: 32idcSP1:)  
--beamline {None,2-bm,7-bm,32-id}    Customized the goodle slide to the beamline selected (default:  
→32-id)  
--idx IDX              Id of x slice for reconstruction visualization (default: -1)  
--idy IDY              Id of y slice for reconstruction visualization (default: -1)  
--idz IDZ              Id of z slice for reconstruction visualization (default: -1)  
--max MAX              Maximum threshold value for reconstruction visualization  
→(default: 0.0)          Minimum threshold value for reconstruction visualization  
→(default: 0.0)  
--presentation-url PRESENTATION_URL          Google presentation url (default: None)  
--rec-type {recgpu,rec}          Specify the prefix of the recon folder (default: recgpu)  
--config FILE          File name of configuration file (default: /home/beams/FAST/logs/  
→tomolog.conf)  
--config-update        When set, the content of the config file is updated using the  
→current params values (default: False)  
--double-fov          Set to true for 0-360 data sets (default: False)  
--logs-home FILE       Log file directory (default: /home/beams/FAST/logs)
```

(continues on next page)

(continued from previous page)

--token-home FILE	Token file directory (default: /home/beams/FAST/tokens)
--verbose	Verbose output (default: False)

For other options:

```
$ tomolog -h
usage: tomolog [-h] [--config FILE]  ...

optional arguments:
-h, --help      show this help message and exit
--config FILE  File name of configuration file

Commands:

init          Create configuration file
run           Run data logging to google slides
status        Show the tomolog status
```

## API REFERENCE

**tomolog-cli Modules:**

### 3.1 tomolog\_cli.tomolog

**Functions:**

---

<i>TomoLog(args)</i>	Class to publish experiment meta data, tomography projection and reconstruction on a google slide document.
----------------------	-------------------------------------------------------------------------------------------------------------

---

```
class tomolog_cli.tomolog.TomoLog(args)
    Bases: object
    Class to publish experiment meta data, tomography projection and reconstruction on a google slide document.

    init_slide()
    plot_projection(proj, fname)
    plot_recon(recon, fname)
    publish_descr(presentation_id, page_id)
    publish_proj(presentation_id, page_id, proj, resolution=1)
    publish_recon(presentation_id, page_id, recon, resolution=1)
    read_meta_item(template)
    read_raw()
    read_recon()
    run_log()
```

## 3.2 tomolog\_cli.tomolog\_2bm

### Functions:

---

<code>TomoLog2BM(args)</code>	Class to publish experiment meta data, tomography projection and reconstruction on a google slide document.
-------------------------------	-------------------------------------------------------------------------------------------------------------

---

```
class tomolog_cli.tomolog_2bm.TomoLog2BM(args)
    Bases: TomoLog

    Class to publish experiment meta data, tomography projection and reconstruction on a google slide document.

    plot_projection(proj, fname)
    plot_recon(recon, fname)
    publish_descr(presentation_id, page_id)
    publish_proj(presentation_id, page_id, proj)
    publish_recon(presentation_id, page_id, recon)
    read_raw()
    read_rec_line()
    read_recon()
    read_tiff_part(fname, x, y, z_start, z0_start, lchunk)
    run_log()
```

## 3.3 tomolog\_cli.tomolog\_32id

### Functions:

---

<code>TomoLog32ID(args)</code>	Class to publish experiment meta data, tomography projection and reconstruction on a google slide document.
--------------------------------	-------------------------------------------------------------------------------------------------------------

---

```
class tomolog_cli.tomolog_32id.TomoLog32ID(args)
    Bases: TomoLog

    Class to publish experiment meta data, tomography projection and reconstruction on a google slide document.

    plot_projection(proj, fname, scalebar='nano')
    plot_recon(recon, fname)
    publish_descr(presentation_id, page_id)
    publish_proj(presentation_id, page_id, proj)
    publish_recon(presentation_id, page_id, recon)
```

```
read_raw()  
read_rec_line()  
read_recon()  
run_log()
```



---

**CHAPTER  
FOUR**

---

**CREDITS**

## 4.1 Citations

GPU based tomographic reconstruction is available at [tomocup](#)

Cite [A1] if you use [tomobank](#)

Cite [A2] if you use [tomopy/tomopy cli](#)

Cite [A3] if you use [DataExchange](#)

## 4.2 References



## BIBLIOGRAPHY

- [A1] Francesco De Carlo, Doga Gursoy, Daniel Jackson Ching, Kees Joost Batenburg, Wolfgang Ludwig, Lucia Mancini, Federica Marone, Rajmund Mokso, Daniel M. Pelt, Jan Sijbers, and Mark Rivers. Tomobank: a tomographic data repository for computational x-ray science. *Measurement Science and Technology*, 2017. URL: <https://doi.org/10.1088/1361-6501/aa9c19>.
- [A2] Gürsoy D, De Carlo F, Xiao X, and Jacobsen C. Tomopy: a framework for the analysis of synchrotron tomographic data. *Journal of Synchrotron Radiation*, 21(5):1188–1193, 2014.
- [A3] De Carlo F, Gursoy D, Marone F, Rivers M, Parkinson YD, Khan F, Schwarz N, Vine DJ, Vogt S, Gleber SC, Narayanan S, Newville M, Lanzilotti T, Sun Y, Hong YP, and Jacobsen C. Scientific data exchange: a schema for hdf5-based storage of raw and analyzed data. *Journal of Synchrotron Radiation*, 21(6):1224–1230, 2014.
- [B1] Nghia T. Vo, Robert C. Atwood, Michael Drakopoulos, and Thomas Connolley. Data processing methods and data acquisition for samples larger than the field of view in parallel-beam tomography. *Opt. Express*, 29(12):17849–17874, Jun 2021. URL: <http://www.opticsexpress.org/abstract.cfm?URI=oe-29-12-17849>, doi:10.1364/OE.418448.



## PYTHON MODULE INDEX

t

tomolog\_cli.tomolog, [17](#)  
tomolog\_cli.tomolog\_2bm, [18](#)  
tomolog\_cli.tomolog\_32id, [18](#)



# INDEX

|

init\_slide() (*tomolog\_cli.tomolog.TomoLog* method),  
17

M

module  
  *tomolog\_cli.tomolog*, 17  
  *tomolog\_cli.tomolog\_2bm*, 18  
  *tomolog\_cli.tomolog\_32id*, 18

P

plot\_projection() (*tomolog\_cli.tomolog.TomoLog* method), 17  
plot\_projection() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
plot\_projection() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 18  
plot\_recon() (*tomolog\_cli.tomolog.TomoLog* method), 17  
plot\_recon() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
plot\_recon() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 18  
publish\_descr() (*tomolog\_cli.tomolog.TomoLog* method), 17  
publish\_descr() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
publish\_descr() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 18  
publish\_proj() (*tomolog\_cli.tomolog.TomoLog* method), 17  
publish\_proj() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
publish\_proj() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 18  
publish\_recon() (*tomolog\_cli.tomolog.TomoLog* method), 17

publish\_recon() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
publish\_recon() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 18

R

read\_meta\_item() (*tomolog\_cli.tomolog.TomoLog* method), 17  
read\_raw() (*tomolog\_cli.tomolog.TomoLog* method), 17  
read\_raw() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
read\_raw() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 18  
read\_rec\_line() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
read\_rec\_line() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 19  
read\_recon() (*tomolog\_cli.tomolog.TomoLog* method), 17  
read\_recon() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
read\_recon() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 19  
read\_tiff\_part() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18

run\_log() (*tomolog\_cli.tomolog.TomoLog* method), 17  
run\_log() (*tomolog\_cli.tomolog\_2bm.TomoLog2BM* method), 18  
run\_log() (*tomolog\_cli.tomolog\_32id.TomoLog32ID* method), 19

T

*TomoLog* (class in *tomolog\_cli.tomolog*), 17  
*TomoLog2BM* (class in *tomolog\_cli.tomolog\_2bm*), 18  
*TomoLog32ID* (class in *tomolog\_cli.tomolog\_32id*), 18  
*tomolog\_cli.tomolog* module, 17

```
tomolog_cli.tomolog_2bm
    module, 18
tomolog_cli.tomolog_32id
    module, 18
```